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Wash 0.1 mmol resin with DMF.

Dissolve 0.244 g (+)-biotin (1 mmol, MW 244.3) in 5 mL DMF-DMSO (1:1) solution. A little warming is necessary.

Add 2.1 mL 0.45 M HBTU/HOBt solution and 0.3 mL DIEA to the solution prepared in step 2.

Add the activated biotin solution to the resin and let stir overnight.

Check resin to make sure coupling is complete as evidenced by negative ninhydrin test (colorless).

Wash resin with DMF-DMSO (1:1) (2x) to remove excess (+)-biotin.

Wash resin with DMF (2x) and DCM (2x).

Let the resin dry before proceeding to cleavage.

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Weigh 10 g 2-chlorotriethyl chloride resin (15 mmol) in a reaction vessel, wash with DMF (2x), swell the resin in 50 mL DMF for 10 min, drain vessel.

Weigh 10 mmol Fmoc-amino acid in a test tube, dissolve Fmoc-amino acid in 40 mL DMF, transfer the solution into the reaction vessel above, add 8.7 mL DIEA (50 mmol), swirl mixture for

30 min at room temperature.

Add 5 mL methanol into the reaction vessel and swirl for 5 min.

Drain and wash with DMF (5x).

Check substitution.

Add 50 mL 20% piperidine to remove the Fmoc group. Swirl mixture for 30 min.

Wash with DMF (5x), DCM (2x), put resin on tissue paper over a foam pad and let dry at room temperature overnight under the hood. Cover the resin with another piece of tissue paper, press lightly to break aggregates.

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Weigh duplicate samples of 5 to 10 mg loaded resin in an eppendorf tube, add 1.00 mL 20% piperidine/DMF, shake for 20 min, centrifuge down the resin.

Transfer 100 μ L of the above solution into a tube containing 10 mL DMF, mix well.

Pipette 2 mL DMF into each of the two cells (reference cell and sample cell), set spectrophotometer to zero. Empty the sample cell, transfer 2 mL of the solution from step 2 into the sample cell, check absorbance.

$$\text{Subs} = 101(A)/7.8(w)$$

A = absorbance

w = mg of resin

Check absorbance three times at 301 nm, calculate average substitution.

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Wash resin with DMF (4x) and then drain completely.

Add approximately 10 mL 20% piperidine/DMF to resin. Shake for one min and drain.

Add another 10 mL 20% piperidine/DMF. Shake for 30 min.

Drain reaction vessel and wash resin with DMF (4x). Make sure there is no piperidine remaining.

Check beads using ninhydrin test, beads should be blue.

Coupling Step - Prepare the following solution:

1 mmol Fmoc-amino acid

2.1 mL 0.45 M HBTU/HOBT (1mmol)

348 μ L DIEA (2 mmol)

Add above solution to the resin and shake for a minimum of 30 min. This coupling step can be longer if desired.

Drain reaction vessel and wash resin with DMF (4x).

Perform Ninhydrin test:

If negative (colorless), proceed to step 2 and continue synthesis.

If positive (blue), return to step 5 and re-couple the same Fmoc-amino acid. Increase the coupling time if necessary.

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For ABI synthesizers, pack Fmoc-Tyr(PO₃H₂)-OH in a cartridge.

The cycle program for coupling Fmoc-Tyr (PO₃H₂)-OH is the same as for other Fmoc-amino

acids except for the coupling time (see step 3). (Note: ABI synthesizers use HBTU/HOBT as the activating reagent.)

The coupling time for Fmoc-Tyr(PO₃H₂)-OH needs to be increased. For ABI model 430A peptide synthesizer, insert several steps (i.e., vortex on, wait 990 sec, vortex off, to increase the coupling time). For ABI model 431A peptide synthesizer, add additional "I"s. Overnight coupling may be necessary for some sequences.

After the coupling step for Fmoc-Tyr(PO₃H₂)-OH, perform ninhydrin test to ensure complete coupling. Negative (colorless) ninhydrin test indicates complete coupling, while a positive (blue) ninhydrin test indicates incomplete coupling.

Increase the coupling time of the amino acid residues after the phosphotyrosine or perform double coupling. (Note: The coupling of amino acids after the phosphotyrosine can be difficult.)

There is a limit on the number of amino acid residues that can be coupled after the phosphotyrosine. Since the phospho group is unprotected, side reactions are likely to occur. (Note: Peptides have been successfully coupled with sequences containing up to ten additional amino acids following the phosphotyrosine residue.)

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Peptides which differ in the C-termini can be simultaneously synthesized in one reaction vessel by employing resins that possess different cleavage properties. The resins used were the weak acid labile 2-chlorotrityl resins and the TFA labile Wang resins. The success of this approach was shown by the co-synthesis of ACTH (4-10) with ACTH (4-11) and Neuropeptide Y, a C-terminal amide peptide with its corresponding C-terminal free acid analog.

*Hong A., Le T., and Phan T. Techniques in Protein Chemistry VI, 531-562 (1995).

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Starting Resin: Chlorotrityl resins

Reagents for 1 g Peptide-Resin:

1 mL acetic acid (AcOH)

2 mL trifluoroethanol (TFE)

7 mL dichloromethane (DCM)

Prepare above mixture.

Add peptide-resin to the mixture and let it stir at room temperature for 1 h.

Filter and wash resin with 10 mL TFE:DCM (2:8) (2x) to ensure that all of the product is recovered.

Evaporate the solvent until there is less than 5 mL of liquid.

Add ether to a test tube containing about 100 mL of the above solution. Check solubility of the fully protected peptide in ether. If the product precipitates, proceed to step 6. If no precipitate is observed, proceed to step 7.

Add cold ether to the residual liquid in step 4 to precipitate the fully protected peptide. Filter through a fine sintered funnel to obtain the product.

Some fully protected peptides are soluble in ether. In this case, add water to precipitate them out. Filter through a fine sintered funnel to obtain the product.

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Fmoc-Lys(Mtt)-OH

Swell resin in DCM.

Wash resin with 3% TFA/DCM (2x) (since the resin is swollen in DCM, this step of washing the resin quickly with 3% TFA/DCM ensures that the actual concentration of TFA is 3%).

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FITC

Fmoc-e-Ahx-OH

Couple Fmoc-e-Ahx-OH to the amino terminal of the peptide-resin using standard coupling conditions.

"De-Fmoc" with piperidine using the standard 20% piperidine procedure.

Wash resin with DMF (3-4x).

Swell resin with DCM and drain.

Prepare solution of 1.1 equivalent of FITC in pyridine/DMF/DCM (12:7:5). Use just enough solution to form a slurry with the resin. Do not use too much solution since the rate of the reaction is proportionate to the concentration of the solution.

Add the solution prepared in step 2 to the resin.

Let mix overnight.

Check the completion of the reaction using ninhydrin test.

If the coupling of FITC to the amino group is not complete, ninhydrin test will give a blue color.

Repeat the coupling with FITC (steps 5-7) if necessary.

Wash resin with DMF (2x), isopropanol (2x), and DCM (2x).

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Reagent:

5-carboxyfluorescein (5-FAM)

Use standard coupling method to couple 5-carboxyfluorescein to the amino group of the peptide. For cost saving purposes, use 2x excess compared to the mmol of resin, instead of the standard 4x excess used for Fmoc-amino acids. For 0.1 mmol synthesis, use 75 mg 5-carboxyfluorescein, 76 mg HBTU, and 70 mL DIEA.